

Compact, reliable and cost-effective

AU480 Clinical Chemistry System



Blood Banking
Centrifugation
Chemistry
Flow Cytometry
Hematology
Hemostasis
Immunoassay
Information Systems
Lab Automation
Molecular Diagnostics
Rapid Diagnostics

The AU480 is the ideal primary chemistry analyzer for small- to medium-sized hospitals and laboratories, or as a dedicated specialty chemical or STAT analyzer for larger institutions. With random access throughput of up to 400 photometric tests per hour (up to 800 with electrolytes), increased onboard testing, reduced sample volume and user-friendly operation, the AU480 delivers efficiency for laboratories around the world.

- New and intuitive graphical user interface, standardized with AU680
 - Sample tracking
 - Metering data
 - Patient statistics
 - User customized menu
 - Color alerts to highlight system operating conditions
- AU proven reliability with quick and easy maintenance procedures
 - No tools required
 - Daily hands on time only 5 minutes
 - Online maintenance videos
- Cooled STAT module provides one button STAT interrupt and advanced Auto QC and calibration capabilities
- High-quality, permanent glass cuvettes—reduces disposable costs
- High-precision micro sampling, ideal for pediatric, geriatric and veterinary testing
- Economical ISEs with long onboard stability; easy to maintain (individual electrode replacement only required)
- 80-sample continuous rack loader